

# Kröber 4.0

- The top model of premium equipment oxygen concentrator Kröber 4.0

## Oxygen concentrator

- super quiet, sensational 31 dB (A) sound pressure level
- reduced weight, only 16 kg
- reduced energy consumption, only 280 Watt



- easy to maintain



- innovative details



- coarse dust filter

# Kröber 4.0

## Oxygen concentrator

- modular components set up
- easily accessible filter and fuses
- USB – interface
- modern design and simple operation
- new, innovative measuring sensor
- integrated flow adjustment
- 30.000 hours warranty for all functional parts
- whisper silent
- microprocessor-controlled
- favourable price
- Made in Germany

### Accessories

- instructions for use
- angular connector for humidifier
- humidifier, refillable
- nasal cannula 2m and 5m

### Technical data

MDD classification	II a
Operating power	230 V 50 Hz, optional 115 V 60 Hz
Temperature	operation +5° to +40°C storage -20° to +50°C
Sound level	31 dB (A)
Power consumption	280 W
Dust filter	in rear of device
Fine filter	behind the maintenance lid
Fuses	mains 230 V 2xTT2.5 A H 250 V optional 115 V 2xTT5.0 A H 250 V internal 1xT1.0 A L 250 V
Alarms/Status indicators	temperature, power failure, flow, oxygen, system
Weight	ca. 16 kg
Dimensions	535x203x520 mm without rollers
Manufacturer's warranty	30.000 operating hours (max. 5 years) according to terms of warranty
Flow adjustment	0 to 2 l/min in 0.1 steps 2 to 4 l/min in 0.2 steps 4 to 5 l/min in 0.5 steps
O2-concentration	1 up to 4 l/min 95% -3% 4 up to 5 l/min 85% -3%



E-mail: [support@baro2.com](mailto:support@baro2.com) Phone: 0345 100 0084

Bar O2 Limited T/A the OxygenStore, Unit 5A Mount Industrial Estate, Mount Road, Stone, Staffordshire, ST15 8LL

Company Registered in UK No 04192727 VAT Registration No 770 4717 24  
EORI GB770471724000 Dun and Bradstreet D-U-N-s No 232439047