



Oxygen simplified

EverFlo oxygen concentrator

The smaller, lighter EverFlo oxygen concentrator from Philips Respironics is making a big impression on providers and patients. Its sleek, ergonomic design makes EverFlo easier to store, move, and carry than conventional concentrators. And it is one of the quietest and most energy efficient devices available today. With no filter for patients to maintain, a durable metal oxygen outlet, and fewer moving parts than other concentrators, EverFlo is designed to deliver consistent, trouble-free operation.





Molded handle that's easy to pick up and transport with one or both hands



Humidifier bottle holder compatible with all bottle styles



Internal tubing storage keeps tubing conveniently available for use



Metal cannula connector for greater durability



Recessed flow meter reduces the risk of breakage



Interface panel with on/off switch and alarm indicators

Product specifications	
	Low power version
Power	230V/50Hz
Liter flow	0.5-5 liters per min
Sound level	43 dBA < 40 dBA
Oxygen purity	87% - 96%
Weight	14 kg (32 lbs)
Dimensions (HxWxD)	580mm (22.8"), 380mm (15"), 240mm (9.5")
OPI (oxygen percentage indicator) alarm levels	Low oxygen: 82%, very low oxygen: < 70%
Required maintenance	Compressor intake filter: 2yrs. Micro disk filter: as needed
Power consumption	< 300 W
Part numbers	1020008 UK, Ireland



PHILIPS
RESPIRONICS

E-mail: support@baro2.com

Phone: 0345 100 0084

Bar O2 Limited T/A the OxygenStore, Unit 5A Mount Industrial Estate, Mount Road, Stone, Staffordshire, ST15 8LL

Company Registered in UK No 04192727 VAT Registration No 770 4717 24
EORI GB770471724000 Dun and Bradstreet D-U-N-s No 232439047